ELEC engine

IMPLANT SURGICAL ENGINE

The ELEC ENGINE is a compact, stylish and equipped with powerful optic motor, real time speed and torque performance exhibited on LCD.







Max 80N.cm

(32:1 gear handpiece)



Torque limit function

The beep sounds when load torque above torque value set by a user operates longer than 3 seconds and the operation stops to ensure patient's safety.



Auto calibration function

Calibration function operates by simultaneously pressing $\blacktriangle/\blacktriangledown$ button for $2\sim3$ seconds.

Perform the Calibration function when the motor is initially used and operate this function under no load status.



Powerful optic motor

It shows LED functioning within motor during the motor operation,







Real-time torque and RPM display

The compact control unit shows real time rpm and torque through digital panel.

It helps users to perform a surgery more precisely.



Program modes

Program button enables users to select 6 modes in turns Drilling, Tapping, Remove Tap, Implant, Remove, Lock Screw. Users can select them for each purpose,



Memory function

Users can save 9 personal settings using memory button, which makes it possible to load the different settings immediately in terms of program mode, rpm, torque, coolant, speed, gear ratio, turn direction.








Sterilizable and Reusable Irrigation tube

Multi-use tubing that can be reused for at least 30 times (with rust-proof cutting edge material)



Thread Cutting Function

Thread cutting is available in Tapping mode among 6 basic programs.



Thread Cutting?

The thread cutter function supports the surgery when makes tap for implants in to dense (D1) bone.



SELECT TAPPING & SET TORQUE



WHEN REACH PRESET TORQUE



IT AUTOMATICALLY REVERSES







MOTOR (ISM -B70L)
Model ISM-B70L
Max.Speed 40,000 RPM
Max.Torque 80 N.cm (32:1)
Internal Resistance 730mQ
Internal Inductance 100uH
Dimension 24mm(W) x 100mm(W)
Weight 156g
Optic Yes

CONTROLLER (ISE-270M)
Model ISE-270M
Power Supply voltage $AC110 \sim 240 \text{V}$
Frequency 50/60Hz
Power Consumption 150VA
Max. Coolant flow rate 110ml/min
Fuse 2 x 250V - T2.0AH
Dimension 285mm(W) x 250mm(W) x 120mm(H)

FOOT PEDAL (FC-V2)
Model FC-V2
Speed control Variable
Foot control Function
Program / Coolant / Rotation Direction







Powerful torque up to 70 N.cm

Reliable power under any operation conditions.



LED brightness up to 25,000 LUX

Excellent light effect can minimize shady spots.



Easy to assemble and disassemble the device

Easy to disassemble and assemble with one touch. handpiece is easy to clean and manage and extend the life of the product



Anti-scratch

Functional design and anti-scratch surface treatment.









Model SG200L
Ratio 20:1
Max.Speed 2,000rpm
Max.Torque 70 N.cm
Size 20mm(D) x 95mm(L)
Weight 73g
Material Stainless steel
Irrigation External
Coolant volume ≥50ml/min
Optic Fiber optic
Warranty 2 Years
Bur CA bur(Ø2.35)
Chuck type Push button
Coupling ISO3963 Type 4
Rotary Instrument IOS1797 Type 1



Model SG200
Ratio 20:1
Max.Speed 2,000rpm
Max.Torque 70 N.cm
Size 20mm(D) x 95mm(L)
Weight 73g
Material Stainless steel
Irrigation External
Coolant volume ≥50ml/min
Optic No
Warranty 2 Years
Bur CA bur(Ø2.35)
Chuck type Push button
Coupling ISO3963 Type 4
Rotary Instrument IOS1797 Type 1



1639

The EU directive 93/42/EEC was applied in the design and production of this medical device.





MICRO-NX Co., Ltd.

22, Maeyeo-ro 1-gil, 41059 Republic of Korea

Tel: 82 53 650 1000 Fax: 82 53 650 1001

MICRO-NX Co., Ltd. branch office

Karl-Marx-Str. 6 16540 Hohen Neuendorf Germany

Tel: 49 (0)3303 5412323 Fax; 49 (0)3303 5412324





Official Website

Official Yutube Chanel