


The soft monofilament

- ▶ Non-wicking monofilament construction
- ▶ Soft and comfortable for patients
- ▶ 300 Series Stainless Steel Needles




Available Sizes


Undyed 19 mm	Precision RC	2/0 USP	CS0418
Undyed 16 mm	Precision RC	3/0 USP	CS0518
Undyed 19 mm	Precision RC	3/0 USP	CS051819
Undyed 16 mm	RC black needle	3/0 USP	CS0518BK
Undyed 19 mm	RC black needle	3/0 USP	CS051819BK
Undyed 13 mm	TP	4/0 USP	CS0618PERIO
Undyed 13 mm	Precision RC	4/0 USP	CS0618PREM
Undyed 16 mm	Precision RC	4/0 USP	CS0618RC
Undyed 13 mm	Precision RC	5/0 USP	CS071813
Undyed 16 mm	Precision RC	5/0 USP	CS071816



Smooth Monofilament Rod



Cytoplast™ PTFE Suture
350x Magnification



PTFE Suture Competitor
350x Magnification

